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Automated External Defibrillator
with Manual Function/ Automated
External Defibrillator

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This quick guide does not substitute the careful reading of the user manual. For complete information about the use of this device, including safety warnings, charging the battery and the correct measurement for the monitored parameters, see the user manual found on the CD that accompanies the equipment.

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Warnings

IMPORTANT: this equipment may only be operated by qualified technical personnel. Read this guide carefully before using the equipment.

WARNING: not recommended for patients younger than one (01) year old.

ATTENTION: the patient's movements can confuse the correct detection of rhythm and delay therapy. Do not perform maneuvers with the patient and keep him still during rhythm analysis.

WARNING: L.on/I.on PRO can be used by patients over 01 (one) year old, regardless of their weight.

WARNING: the patient must be placed on non-conductive surfaces. Do not use wet or metallic surfaces and, if necessary, dry the chest before applying the shock.

WARNING: do not touch the patient, the equipment, the accessories nor any metallic or conductive surface which is in contact with the patient during the defibrillation.

WARNING: the patient must be completely still during the cardiac rhythm analysis phase. Do not give cardiac massage at this point.

WARNING: risk of explosion if the equipment is operated in the presence of flammable liquids or gases.

WARNING: always check the general state of the equipment, the battery and the accessories before using it.

NOTICE: each and every repair to the equipment can only be done by intramed's authorized technical assistance centers.

NOTICE: the use of the L.on/I.on PRO is restricted to one patient at a time.

NOTICE: the applied parts are protected against defibrillation discharge.

NOTICE: during discharge there may be baseline variation.

NOTICE: the limits of current leakage may be exceeded.

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NOTICE: avoid connecting the patient to several items of equipment at the same time.

NOTICE: the applied parts intended to come into contact with the patient have been evaluated and comply with the directives and principles of ISO 10993-1.

NOTICE: when removing the equipment from the package, carefully verify if there is any abnormality or visible damage in the device or its accessories, caused by impact or mishandling during transportation. In case of irregularities, please contact Instramed.

NOTICE: disposable accessories and any other components must be disposed of according to the norms of hospital waste disposal.

NOTICE: do not modify this equipment without authorization from Instramed.

ATTENTION: the equipment may not be connected to the patient when communicating via USB with the SoftIDEA app.

ATTENTION: the equipment blocks patient operation when communicating via USB with a PC.

WARNING: the I.on/I.on PRO is designed to have no loose parts or small parts, but if any parts come loose, there is a risk of asphyxiation from ingestion or inhalation. Therefore, keep the equipment and its materials out of the reach of children.

WARNING: The I.on/I.on PRO and its accessories are free of latex and allergy-causing components.

ATTENTION: the I.on/I.on PRO may suffer interference or its performance may be affected if it is near sources of heat or moisture, such as near heaters, kitchen equipment, or open areas. Therefore, try to keep I.on/I.on PRO in areas protected from these conditions to ensure its perfect operation.

Attention

The following factors may cause ECG misinterpretation:

- Wrongly placed pads.
- Patient's movements.

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- Pacemaker (it may lessen the precision of the cardiac arrest detector).
- Radio frequency interference, including mobile phones.
- Excessive hair or wet skin in the application area of the electrodes.
- Pieces of clothing between skin and pads.

I.on/I.on PRO ONLY works with battery.

Warnings, precautions and restrictions

The I.on/I.on PRO is a device built according to NBR and IEC standards and therefore is completely safe for the patient and the user. However, all safety precautions described below must be followed.

The operation of the I.on/I.on PRO may be affected by the presence of electromagnetic power supplies, such as electrosurgical equipment and computer tomography (CT).

Electromagnetic compatibility

Warnings

Using the I.on/I.on PRO requires special precautions concerning Electromagnetic Compatibility in compliance with the information contained in this manual.

Mobile and portable RF communications equipment, such as a cell phone, may affect the functioning of the I.on/I.on PRO.

Maximum length of accessory cables in order to comply with the requirements of Electromagnetic Compatibility is 2.5 m.

All parts and accessories, listed below, follow the requirements for electromagnetic compatibility.

- Pair of disposable adhesive pads, adult size.
- Pair of disposable adhesive pads, child size (optional).
- Power supply for charging the battery (optional).
- USB cable.

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Warnings

Using cables and accessories different from the ones specified above, except for cables and accessories sold by Instramed as replacement pieces, may result in emission gain or immunity decrease of the equipment.

The I.on/I.on PRO must not be used too close to or piled over other equipment.

The I.on PRO is suitable for use in professional health care environments.

Do not use the equipment in the presence of MRI devices. Measures must be taken to minimize interference from high frequency surgical equipment.

Portable RF communication equipment (including peripherals such as cables and external antennas) should not be used within 30 cm of any part of the I.on PRO, including cables specified by Instramed. Otherwise, performance degradation of this equipment may occur.

Adverse or side effects

Superficial burns may occur to the patient's skin on the area in contact with the electrodes. To minimize this effect, apply the pads soon after removing their protective envelope and press them firmly to the patient's skin.

Possibility of superficial skin burns. To minimize the effect, in the case of adhesive pads, apply them immediately after removal of the protective envelope and securely attach to the patient's skin. The patient's skin must be dry.

Possibility of reduction of treatment efficiency. The patient's skin must be dry, otherwise the electric discharge may leak. Never apply conductive gel.

Possibility of reddish and/or bruised skin at the application place (thorax) by the use of CPR MAESTRO. It is recommended for cases of resuscitation maneuvers of long duration, the use of a gauze between the skin and the CPR MAESTRO.

Possible BURNED/REDNESS SKIN, due to the HIGH VOLTAGE and HIGH CURRENT delivered. But it becomes more severe damaged as the progressive counting of the delivered shocks. These side-effects depend also on the quantity of delivered shocks.

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ECZEMAS on the skin, due to a bad biocompatibility of the adhesive pads or ECG electrodes. The accessories supplied with the product are biocompatible according to the ISO 10993 standard.

Degraded or loose electrodes may affect equipment performance, possibly resulting in reduced efficiency or inability to treat.

The effects of animals or pests, such as birds, rodents, among others, can cause perforation, broken connections or loss of contact, impairing the performance of the equipment. Therefore, try to keep l.on/lon PRO in areas protected from these conditions to ensure its perfect operation.

The l.on/lon PRO has bright colors to facilitate identification in situations of need. On the other hand, it may happen to attract the attention of children. To avoid the risk of choking or degrading the performance of the equipment, keep it out of the reach of children.

Device care

- Do not place the equipment where it may fall on the patient. Do not lift the equipment by its cables or connections.
- Place cables connected to the patient in order to restrict the possibility of strangulation.
- Always keep the equipment and its accessories clean and well maintained.
- If you suspect a fall or external damage, do not use the equipment.

Important recommendations

- Never immerse in liquid and never spill liquid of any kind on any part of the equipment.
- Do not use any other cleaning products not recommended by this manual.

ATTENTION: NEVER sterilize any parts of the equipment, regardless of the sterilization method, as this would damage the mechanical structure and compromise the product's operation.

The above recommendations ensure that the device will support, without damage or deterioration of safety factors, the cleaning and disinfection process necessary.

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Qualified users

Qualified users will be considered those who have been trained in a recognized training school in the use of automatic defibrillators and CPR - Cardiopulmonary Resuscitation techniques. Although the new regulations of July 2021 do not provide for specific training under the "Good Samaritan" law.

Disposing of the device

As provided in Brazilian environmental legislation, equipment and parts that do not have more use conditions should be referred to the manufacturer for proper final destination, thus preserving natural resources and contributing to the conservation of the environment.

For disposal of Instramed brand products, contact us through the phone numbers available at www.medicalfast.eu

Avoid contamination of the environment, humans, or other equipment by make sure to properly decontaminate the equipment before disposing of it.

Refer to local regulations for the proper disposal of trash in your area. For countries that follow European Guidelines, refer to 2002/96/CE.

More information on this topic is available in the Ion/Ion PRO User Manual present on the CD which accompanies the product.

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Included items

- An I.on/I.on PRO automated defibrillator.
- A pair of disposable adult size adhesive pads.
- One first-aid kit containing 1 pair of surgical gloves, 1 scissor and 1 CPR mask.
- A transport bag.
- A USB cable.
- CD with instructions manual and SoftDAE management PC software.

Optional items

- A power supply for charging the battery.
- A cable for connecting the power supply to the electric current.
- A pair of disposable child size adhesive pads.
- 3 lead ECG cable.
- Accessory for chest compression (CPR Maestro)

Replacement parts

Request MedicalFast to replace the following items:

- Batteries replacement.
- Adult/child adhesive pads replacement.

To request pieces and services, contact MedicalFast srl.

 More information on this topic is available in the I.on/I.on PRO User Manual present on the CD which accompanies the product.

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To charge the battery, disconnect it from the I.on/I.on PRO and connect the charger at the back of the battery and then connect it to the mains supply.

- Time to battery full charge = 5 hours.

ATTENTION: the equipment may not be connected to the patient when communicating via USB with the SoftDEA app.

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ATTENTION: I.on/I.on PRO devices only operate on battery power.

Occasional use

Even when disconnected (stand-by), the I.on/I.on PRO executes internal routines checking the status of the equipment. Despite of this procedure entailing a low power consumption, the battery charge may be consumed. Therefore, it is recommended that the battery be fully charged every 8 months.

Rechargeable batteries replacement

Every battery has a determined shelf life, which is the possible quantity of full charge and discharge cycles, without loss of performance (see battery specifications in chapter 8). When the device presents a drop in battery performance, with low autonomy, please request a new unit from Instramed's technical assistance.

The battery can be withdrawn through the rear opening through a quick-release system. Remove the old battery and replace the new battery, observing the correct docking position.

Battery replacement is recommended every 02 (two) years or when the runtime is less than 01 (one) hour.

More information on this topic is available in the I.on/I.on PRO User Manual present on the CD which accompanies the product.

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Front panel


I.on and I.on PRO models

- LCD screen version for health care professional.

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1. Touch screen: presents operational information and allows for manual interaction with the device.
2. Operational status indicator.
3. Speaker.
4. Start button.
5. Microphone (optional).
6. Patient selection button for health care workers (optional).
7. Handle set (optional).

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
Modello I.on

- Con schermo LCD.

Touch screen for health care professionals

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
1. Heart beats per minute.
2. CPR interval counter: counts the interval between discharge delivery, helping in the CPR.
3. Shocks counter: shows the total number of defibrillations successfully executed.
4. Battery status.
5. Orientation message.
6. ECG curve.
7. Icon indicating the defibrillation stage.

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1. Heart beats per minute.
2. Shocks counter: shows the total number of defibrillations successfully executed.
3. General timer: shows the total time of the equipment being on.
4. CPR interval counter: counts the interval between discharge delivery, helping in the CPR (Cardiopulmonary resuscitation) massage.
5. Battery status.
6. Orientation message.
7. Manual mode button.
8. ECG curve.
9. Icon indicating the defibrillation stage.

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
1. Heart beats per minute.
2. Shocks counter: shows the total number of defibrillations successfully executed.
3. General timer: shows the total time of the equipment being on.
4. CPR interval counter: in manual mode this indicator shows the continuous counting since the beginning of the operation.
5. Battery status.
6. Charge selectors: use these buttons to adjust the charge in J.
7. Selected charge.
8. ECG curve.
9. Automatic mode button.
10. Charge button: charges the defibrillator to the selected energy setting.
11. Cancel button: disarm the stored charge.

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
The assistance step is shown by the LED indication above the pictogram.

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Operational status indicator

The I.on/I.on PRO performs a weekly complete auto test as standard* allowing the user to know the operational status of the device. This status is informed by a visual indicator (see pictures below), voice messages and sound signals. The automatic test is also performed when the device is switched on.

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VISUAL INDICATOR

The visual indicator LEDs flash at a 5-second interval to report the following status.

- Indicates that the equipment is operational and ready for use.

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- It indicates that the device **DOES NOT HAVE ENOUGH CHARGE IN THE BATTERY TO OPERATE** or presents another internal defect. In case of non-rechargeable battery, change the battery immediately. In case of rechargeable battery, charge the battery immediately. After changing or recharging the battery, turn on the equipment to perform the auto test. If the indicator remains red, contact Instramed Technical Support or an authorized representative. It indicates that the device **DOES NOT HAVE ENOUGH CHARGE**.

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NOTE: EVEN AFTER THE BATTERY HAS BEEN COMPLETELY CHARGED, the operational status indicator may continue to show an  or some time (rechargeable battery version).

The display will only change from  to  when the l.on/l.on PRO performs the auto test routine or if the device is turned on/off by the user.

 **ATTENTION:** remember to check the status of the operational status indicator at least every 30 days.

SOUND INDICATOR

Along with visual indication, the l.on/l.on PRO emits electronic beeps.

 **ATTENTION:** the device will not turn on in case of low battery or presenting general failure. In this case only the audible beep warning will be issued.

Speaker

The l.on/l.on PRO is a highly complex equipment which, from the moment of activation, assesses the steps of the operation and the general state of the patient. Based on this analysis, the device guides the user through verbal commands which may be warnings, instructions or status messages.

Therefore it is extremely important that the speaker is unobstructed and the l.on/l.on PRO is in a position which allows the user to hear its instructions.

 **ATTENTION:** do not use the equipment inside bags which may prevent the user from hearing the spoken instructions.

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
Start button

The I.on/I.on PRO offers a unique technology that allows the operation of the device to be performed with just one button, completely safe.

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The start button has the functions of:

1. Turning on the device.
2. Starting the automatic process of the patient's clinical analysis.
3. Applying shock therapy (active only when the automatic clinical analysis of the patient indicates the need for it).

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NOTE: it is not necessary to switch the I.on/I.on PRO off. 45 seconds after the removal of pads from the patient or disconnection of pads from the equipment, the device switches itself off, saving battery charge. In this moment the following message will be heard: "The device is being turned off. Press the button in order to turn the machine back on".

Even so, there are two ways to manually switch off: press the start button for three seconds and remove the pads (after 45 seconds without the pads the equipment will automatically shut down).

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Microphone (optional)

The I.on/I.on PRO has the functionality of storing ambient sound. The maximum storage capacity of ambient sound is 10 hours. Stored audios can be transferred to the computer using SoftDEA.

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Patient selection button (optional)

When using electrodes intended for adults, it is possible, through the patient selection button, to switch the equipment's operating mode to adult or child.

Selection of the patient type is indicated by the LEDs in the center of the pictograms, followed by the Adult mode/Child mode confirmation message.

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1. Adult mode.
2. Child mode.

If the I.on/I.on PRO identifies the child electrode connection, the child mode is automatically selected and the patient selection button is inhibited from working.

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1. Disposable pads connector.
2. CPR Maestro connector.
3. ECG connector (only in the models with LCD screen).

ATTENTION: whenever the pads set is replaced, remember to keep the new pair already connected.

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ATTENTION: disposable pads have defined expiration date. Check the enclosure for limit date for use and, if it is not used within this time, replace them with a new pair.

ATTENTION: only use original pads, supplied by Instramed. Failure to do so may result in malfunction.

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
1. Battery compartment.
2. USB connector.

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In case of battery replacement, use original replacements from Instramed supplied by its authorized distributors.

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Charging the battery

To recharge the rechargeable battery, only remove it from the equipment and connect it to the charger, using the input indicated below.

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1. Battery charger connector (rechargeable version)

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
VISUAL INDICATORS

Visual indicators of the charging process may be found in the charger.

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1. Battery LED

- **On**

Indicates that the battery is being charged.

- **Off**

Indicates complete charge.

2. Mains LED

- **On**

Indicates that the equipment is connected to the mains supply.

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
When in AED mode (Automated External Defibrillator), the l.on/i.on PRO identifies arrhythmias and selects the energy charge automatically. l.on/i.on PRO, on AED mode, is in accordance with American Heart Association 2015 Guidelines.

Below you will find a simplified introduction to how the l.on/i.on PRO operates. Be sure to carefully memorize the detailed guide on the next pages before operating the equipment.

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⚠ Before starting the operation, call the emergency service.

1. Press button to start cycle.
2. Place pads on patient.
3. Equipment analyses patient's condition.
4. Analysis is conducted: shock is recommended!
5. Press button to deliver shock.
6. Shock delivered: start CPR. (Cardiopulmonary Resuscitation)

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The energy delivered is pre-adjusted at the factory according to the values below. The operator can only change this protocol using SoftDEA.

- For adult discharge: 1st discharge 200 J - 2nd discharge 300 J - 3rd discharge 360 J successive 360 J.
- For child discharge: 50 J.

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After disposing of used pads, always leave a replacement pair already connected to the equipment, avoiding having to replace them at the moment of the emergency.

ATTENTION: this device has electronic safeguards and will not operate in inadvisable situations.

ATTENTION: check patient's condition. Only use the equipment if the patient is not breathing.

ATTENTION: the area in contact with the pads must be dry.

ATTENTION: the presence of too much hair in the contact area may affect scanning. In this case, shave hair.

ATTENTION: the pads must be applied directly over the skin. DO NOT place pads over clothes.

ATTENTION: the pads are disposable and for single use, cannot be reused under any circumstances.

ATTENTION: after opening the wrapping, the pads should be used within 24 hours.

ATTENTION: in case of use for long periods, the pads should be replaced every 24 hours.

ATTENTION: the patient must be on a steady surface. Any movement during the process of clinical analysis will result in mistaken scans.

ATTENTION: the pads are disposable and can be used in only one patient at a time. Remember to always keep extra ones with the equipment. For replacements, please contact MEDICALFAST.

ATTENTION: to avoid analysis error, no one should touch the patient during analysis.

ATTENTION: the user must not touch the patient or conductive surfaces in contact with him/her during shock delivery, under risk of suffering a powerful electric discharge.

ATTENTION: disconnect other equipment which do not have defibrillation protection before defibrillating the patient.

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Only I.on PRO


⚠ This type of action must be done by a health care professional.

The I.on/I.on PRO enables manual mode operation, as a conventional defibrillator. In this situation, the device does not interfere with the treatment, and the user is responsible for choosing the energy, charging the charge and delivering the shock. After confirming the mode change, I.on/I.on PRO ceases to emit sound and visual orientation, in addition to the automatic safeguards against shocks.

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1. Press button to start cycle.
2. Place pads on patient.
3. Activate manual mode.
4. Select charge and charge energy.
5. Press button to deliver shock.
6. Shock delivered: start CPR. (Cardiopulmonary Resuscitation)

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ATTENTION: the use of the manual mode is the user's entire responsibility. The use by non-qualified professional may cause severe damage and even the patient's death.

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Preventive maintenance

Instramed recommends that the equipment be examined by a qualified technician every 12 months. We recommend that you contact the manufacturer for more information about qualified and trained personnel in your area to perform preventive maintenance.

It is recommended that periodic inspections be performed on the equipment's power supply charger, cables and connectors in order to determine possible isolation or internal conductor ruptures.

Remember to check the status of the operational status indicator at least every 30 days.

Corrective maintenance

If the equipment needs repair, this can only be done by INSTRAMED or its authorized representative, otherwise this Warranty certificate may no longer be valid.

No internal parts are to be fixed by the user.

ATTENTION: periodic maintenance is needed independently of the equipment's use frequency.

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Instramed Indústria Médico Hospitalar Ltda. warrants the equipment described in this Certificate for 12 (twelve) months, starting from the delivery date. This warranty covers manufacturing or material defects that prevents proper functioning according to the specifications stated herein, as long as the conditions presented in this certificate are respected.

During the warranty period, Instramed Indústria Médico Hospitalar Ltda. or its representative will repair or replace defective parts, at no expense to the equipment's owner.

This warranty will no longer be valid if any damage occurs due to accident, natural disaster, improper connection to a power source, use distinct from that described in the User manual, or irregular working conditions.

Any attempt to violate, adjust or repair this equipment by individuals not authorized by Instramed Indústria Médico Hospitalar Ltda. will automatically invalidate this warranty. This also applies in case of alterations made to this contract, the fiscal receipt, or to the equipment's serial number.

Instramed Indústria Médico Hospitalar Ltda. is not responsible for the improper use of this equipment, by people who are not familiar with its function or the techniques recommend for its proper use.

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Automated External Defibrillator with
Manual Function/ Automated External
Defibrillator

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ATTENTION: Instramed assumes no responsibility for any damage caused to individuals or property brought by failure to use this product in accordance with the information, recommendations and warnings presented in the user manual, alterations made in the device, attempts of repair not provided by authorized technical assistance centers, operation by unqualified personnel, use of defective device or use of accessories and parts not supplied by the manufacturer.

For information about warranty or technical assistance, please contact Instramed's technical support.

Copyright © 2019 Instramed. The I.on/I.on PRO, Instramed and its respective logos are trademarks of Instramed Indústria Médico Hospitalar Ltda. The internal software of this product is Instramed's intellectual property, being protected under international copyright laws. It is provided exclusively to be used with this present device, identified by the serial number, and may not be, in whole or in part, evaluated, recompiled or altered in any way.

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