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This quick guide does not substitute the careful reading of the user manual. For complete information about the use of this device, including safety warnings, charging the battery and the correct measurement for the monitored parameters, see the user manual found on the CD that accompanies the equipment.

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English

Attention

-  The following factors can cause bad contact, causing burns on the patient.
- Misplaced pads
- Excessive hair or wet skin in the area of the electrodes application
- Pieces of clothing between skin and pads

Warning

 **CAUTION:** to perform a direct defibrillation (not synchronized), the synchronized defibrillation indicator LED must be off. Otherwise, the defibrillator will not apply the energy discharge on the patient due to absence of ECG signal or electrical connection between the defibrillator and a cardiac monitor for acquisition of the ECG signal.

For synchronized (cardioversion) discharge, the device must detect the ECG signal with QRS.

IMPORTANT: this device must only be operated by qualified technical personnel. Before using it, read the user manual attentively.

ATTENTION: risk of explosion if the equipment is operated in the presence of flammable liquids or gases.

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ELECTRICAL SHOCK HAZARD: never open the device. Each and every repair must be performed by Instramed's authorized technical centers.

ATTENTION: THE PATIENT MUST BE PLACED ON NON-CONDUCTIVE SURFACES. DO NOT USE WET OR METALLIC SURFACES AND, IF NECESSARY, DRY THE PATIENT'S CHEST BEFORE SHOCK DELIVERY.

ATTENTION: DURING THE DEFIBRILLATION, DO NOT TOUCH THE PATIENT, THE EQUIPMENT, THE ACCESSORIES NOR ANY METALLIC OR CONDUCTIVE SURFACE IN CONTACT WITH THE PATIENT.

Do not use Apolus in the presence of MRI equipment.

This equipment was projected to offer resistance to electromagnetic interferences. However, the functioning of this device can be affected in the presence of strong sources of electromagnetic-interference or radio-frequency, such as mobile phones, communicator radios, etc.

ATTENTION: always check the general state of the equipment and its accessories before using it.

Before installing the equipment verify if there are any abnormalities or damage caused by mishandling during transportation.

WARNING: the use of the Apolus is restricted to one patient at a time.

WARNING: the conductive parts of the electrodes and connectors associated with the applied parts, including the neutral electrode, must not come into contact with other conductive parts, including the ground wire.

WARNING: avoid connecting the patient to several equipments at the same time, because the limits of current leakage may be exceeded.

WARNING: in general, the parts of the EQUIPMENT and ACCESSORIES of the Apolus defibrillator intended to come into contact with biological tissues, cells or fluids are tested and analyzed according to the directives and principles of ISO 10993-1, which deals exclusively with the biocompatibility test of the applied parts.

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WARNING: when removing the equipment from its package, carefully verify if there is any abnormality or visible damage in the device or its accessories, caused by impact or mishandling during transportation. In case of irregularities, please contact Instramed.

WARNING: disposable accessories and any other components must be disposed according to the hospital waste disposal norms.

Adverse effects

Superficial burns may occur on the patient's skin in the area in contact with the electrodes. To minimize the effect of the disposable paddles, apply them directly after removal from the protection envelope and attach them firmly to the patient's skin.

The skin must be dry, or electric current leakage may occur, increasing the burn's area and reducing the efficiency of the treatment.

Device care

Do not place the equipment where it may fall on the patient. Do not lift the equipment by its cables or connections.

Place cables connected to the patient in order to restrict the possibility of strangulation.

Keep the defibrillator in a dry environment, avoiding places that allow liquids to spill over the monitor. Do not use the defibrillator if it is wet or excessively humid.

Always keep the equipment and its accessories clean and well maintained.

If you suspect a fall or external damage, do not use the equipment.

Grounding

⚠ GROUNDING IS ESSENTIAL TO PROTECT THE OPERATOR AND PATIENT AGAINST ELECTRICAL DISCHARGE ACCIDENTS. IN THE ABSENCE OF ADEQUATE GROUNDING, DANGEROUS CURRENTS MAY CIRCULATE FROM THE EQUIPMENT BOX IF THERE IS AN INTERNAL ELECTRICAL DEFECT. GROUNDING MUST BE PERFORMED ACCORDING TO ABNT NORMS FOR ELECTRICAL INSTALLATIONS (NBR 13534/1995).

The potential equalization is performed by the Mains Supply cable connector with 3 pins (accompanying the product) and / or with the earthing cable (optional).

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Electromagnetic compatibility

! The installation of the Apolus requires special precautions concerning Electromagnetic Compatibility in compliance with the information contained in this manual.

 More information on this topic is available in the Apolus User Manual present on the CD which accompanies the product.

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- Operational: 0 to 50 °C.
- Storage: 0 to 70 °C.

Humidity

- Operational: 10 to 95% RH, without condensation.
- Storage: 10 to 95% RH, without condensation.

! WARNING: if the Apolus is used outside these conditions, 15 through 30 minutes will be required to stabilize the system so that functioning failures do not occur.

 More information on this topic is available in the Apolus User Manual present on the CD which accompanies the product.

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Included items

When opening the package, please check whether all items below are present:

- Apolus Biphasic Defibrillator.
- Instructions manual.
- 3-pin professional power cable.
- External defibrillation pads set | Adult and child.

Replacement parts

You can call Instramed for replacements of consumable items, parts and accessories. Consult Instramed for prices. Shipping may apply.

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To request pieces and services please contact the representative of your region (the list may be found on www.instramed.com.br) or directly Instramed by the telephone +55 51 3073-8200.

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

1. Power, battery charge and QRS indications.
2. LCD display.
3. Selector switch: turns the equipment on and off; selects the energy.
4. Operation buttons.
5. Transport handle.

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1. **Scale from 5 to 360 Joules:** allows the user to turn the equipment on and select the desired energy.
2. **Position "Off":** turns off the equipment.
3. **Position "Test":** functional test.

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
LCD screen

The LCD screen of the Apolus shows numeric information related to the defibrillation process, as well as the battery level indicator (5 levels).

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Power, battery charge and QRS indications

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1. **Power connected:** when the LED is on, it indicates that the equipment is connected to a power source or an external battery.
2. **Battery charging:** when the LED is on, it indicates that the battery is charging.
3. **Systole visual indication:** this LED lights up during each heartbeat, indicating the recognition of the R wave peak.


NOTE - The LEDs light up even with the equipment is turned off.

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

Operation buttons


Sync: used to enable the synchronized discharge. When the LED is on, it indicates that the function is active. At the same time, the QRS LED blinks according to the heartbeat.

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

Charge: its activation triggers the accumulation of the internal energy that will be used in the shock treatment.

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

Shock: its activation releases, through the pads, the stored energy, resulting in an electric pulse applied to the patient's heart.

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Disarm: disarms the stored internal energy after activation of the "charge" button.

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
1. Connector for defibrillation electrodes (pads)

- **Adult/child external:** included with the equipment, may be used on adults and/or children.
- **Internal:** used in surgeries.

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
1. ECG signal socket: ECG input connector (1V/mV) for synchronized defibrillation.

2. Grounding and potential equalizer: potential equalization and general grounding connector.

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3. **External DC socket:** for battery connection or external DC source in a range of 11 to 16 VDC.
4. **3-pin power connector:** input of 100 to 265 VAC, with central pin for grounding. 5A fuse (20mm 20AG F5A GLASS FUSE).
5. **Battery on/off switch:** internal battery on/off switch. Used to turn off the battery when the equipment is not used for a long period.

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Warnings

- !** The Apolus has a patient impedance meter, it delivers shocks in impedances of 25 to 300 Ohms.
- If a cable or conductor is suspected of being ruptured, avoid using the equipment due to possible risk to the operator.
- Ensure that the defibrillation electrodes of the Apolus are at an appropriate distance from other electrodes so that the power applied does not flow through these electrodes.
- Disconnect all equipment devoid of protection against the discharge of defibrillators.
- Ensure that the patient does not come into contact with any metallic parts.

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

A - Sternum
B - Apex

1. Check if the pads are connected to the Apolus. If they are not, connect the defibrillation cable to the pads socket located on the equipment's side. Turn the screw all the way.
2. Take both pads from their base pulling them up and out.
3. Apply the conductive material to the pads' electrodes.
4. Place the pads on the patient.

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
⚠ ENSURE that the electrodes are away from each other. DO NOT apply paste or gel to the thorax between the pads or the current may follow a superficial route along the thorax wall and not reach the heart.

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⚠ All LEDs blinking indicate short circuit in the pads. Shock delivery will not be allowed.

The red LED blinking indicates bad contact with the patient. Shock delivery will not be allowed.

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The Apolus will automatically identify that it is operating in pediatric mode. Energy is limited to 50 Joules in the pediatric mode.

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Defibrillation

• Step 1 - SELECT ENERGY

Rotate the selection switch until you reach the energy desired. Energy options go from 5 to 360 Joules. In most cases, 200 Joules is recommended for adult use.

• Step 2 - CHARGE

Press the "Charge" button (green in the front panel or use the charge button in the external pads (orange). While the Apolus is charging, a sound will be emitted and the measurement of the charged energy will appear on the display.

The energy selected can be increased or decreased at any time just by rotating the selector switch to the new charge.

To cancel the shock press "DISARM".

When the charge is complete, the device sends a sound signal and displays "Charge Ready" on the screen.

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• Step 3 - SHOCK

After the "Charge Ready" warning, press the "Shock" 3 button (orange) in the front panel or use the two buttons (orange) in the external pads.

Only with the adult/child external pads it is possible to defibrillate using the pad buttons.

CAUTION: make sure nobody is touching the patient! Clearly warn everyone to stand clear of the patient!

The number of shocks and length of operation are indicated on the Apolus display.

Synchronism - Synchronized discharge - Cardioversion

Remember: The function "Synchronized Shock" is disabled after the shock is delivered.

For synchronized defibrillation, a cardiac monitor with ECG output of 1 V/mV is required.

1 - Using a sync cable, connect the ECG input, located on the rear panel of the Apolus, to the 1 V/mV ECG output of your cardiac monitor.

2 - Press the "Sync" button on the panel for two seconds. The green indication next to the button will light up.

3 - Check if the QRS signal on your monitor triggers the Apolus synchronism circuit by checking for a visual (QRS LED) and audible QRS indication on the defibrillator.

In case the Apolus is unable to identify the synchronism pulse, the following message will appear on the display of the device:

NO SYNCRONISM PULSE

Check the connections and cable condition. If this does not resolve, it may be necessary to adjust the gain of the ECG input amplifier on the cardiac monitor. Consult the user manual or the monitor manufacturer.

IMPORTANT: keep key 3 (shock) or the two shock buttons on the pads pressed until the next "R" wave is identified. The Apolus will deliver the shock when the next "R" wave is identified.

IMPORTANT: if Apolus does not identify a valid QRS it will not trigger the shock!

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Disarm key

Disarm the stored charge. Charge may be disarmed at any time, whether the charge is ready or not.

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1. **Energy selected:** from 5 to 360 Joules.
2. **Type of defibrillation electrode:** adult pad, internal pad, adhesive pad or disconnected pad.
3. **Energy charged:** after selection and charging of the energy desired, the following message appears on the display: **CHARGE READY - PA ADULT.**

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Attention


The functional test must be performed daily, which ensures that the equipment is working perfectly and in a state of readiness.

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• Step 1

Put the selector switch in the functional test position.

• Step 2

In case the pads are not connected, the display will show the message: CONNECT THE PADS ON THE EQUIPMENT. Follow the instruction

• Step 3

In case the pads are not positioned in the support, the display will show the message: CONNECT THE PADS ON THE SUPPORT. Follow the instruction.

• Step 4

The display will show the message: PRESS CHARGE. Press the "charge" key and wait until the equipment sends the charge ready signal.

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1. Equipment passed functional test.
2. Equipment failed functional test.

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Attention

 If Apolus fails the functional test, contact technical support immediately.

Preventive maintenance

Instramed recommends that the Apolus be examined by a qualified technician every 12 months. We recommend that you contact the manufacturer for more information about qualified and trained personnel to perform preventive maintenance.

It is recommended that periodic inspections be performed on the equipment's power supply charger, cables and connectors in order to determine possible isolation or internal conductor ruptures.

It is recommended that functional tests be performed at the beginning of every work shift.

Corrective maintenance

Shall be considered qualified users those who have had training in a recognized institution in the use of automated defibrillators and CPR techniques - Cardiopulmonary Resuscitation.

Attention

 No internal parts are to be fixed by the user.

Cleaning and disinfection

Instramed recommends cleaning and disinfecting the equipment and its accessories every three months, or shorter periods whenever excessive dirt or contamination is noticed.

Equipment

- Remove the equipment from the electric current before cleaning it.
- Wipe the external part of the equipment with a cloth dampened with water and neutral soap or isopropyl alcohol.

Accessories

- Use a cloth dampened with isopropyl alcohol.
- Rub the surface to be cleaned for about 10 minutes.

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Sterilization

Instramed recommends that the set of internal pads and their adult electrodes be sterilized using the "cold sterilization" method, making the right use of chemical ETO (ethylene oxide) ensuring quality control of sterilization process and handling by specialized professionals. However, the connector of the internal pad set (which is connected to the equipment) should not be immersed with the entire length of this accessory.

NEVER sterilize any parts of the equipment or its accessories using "dry heat sterilization", such as when using an autoclave. This will damage the mechanical structure and compromise functioning.

Battery

If Apolus is not used for a long period of time, the battery will need to be recharged. To recharge the battery, connect the monitor to an AC power source (110 or 220V outlet) or a DC power source.

There are no restrictions or limitations for using the Apolus while its battery is being recharged by an AC source or DC External source.

Returning components

If Apolus must be returned for repair, call Instramed for shipping instructions. Be prepared to provide the equipment's series number.

If possible, use the original equipment's packaging. If this is not possible, use an equivalent box that provides adequate protection for the monitor.

Precautions, restrictions and warnings

The Apolus is a device built according to NBR and IEC standards and therefore offers total safety for patient and operator. However, all safety precautions described below must be followed.

The equipment may have its operation affected by the presence of electromagnetic power sources, such as electrosurgical equipment and computer tomography (CT).

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Disposing of the device

According to the Brazilian environmental legislation, equipments and parts that are no longer in conditions of use should be referred to the manufacturer for the final destination, thus preserving the natural resources and contributing to the conservation of the environment.

For disposal of products from Instramed, contact us by the telephone numbers available on the website "www.instramed.com.br" or by the e-mail "qualidade@instramed.com.br".

To avoid contamination of the environment, humans, or other equipment, make sure to properly sterilize and decontaminate the equipment before disposing of it.

For countries that follow European Guidelines, refer to 2002/96/CE. For other, countries, refer to local regulations for the proper disposal of trash in your area.

More information on this topic is available in the Apolus User Manual present on the CD which accompanies the product.

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Garantia | Warranty certificate | Garantia

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Instramed Indústria Médico Hospitalar Ltda. warrants the equipment described in this Certificate for 12 (twelve) months, starting from the delivery date. This warranty covers manufacturing or material defects that prevents proper functioning according to the specifications stated herein, as long as the conditions presented in this certificate are respected.

During the warranty period, Instramed Indústria Médico Hospitalar Ltda. or its representative will repair or replace defective parts, at no expense to the equipment's owner.

This warranty will no longer be valid if any damage occurs due to accident, natural disaster, improper connection to a power source, use distinct from that described in the User manual, or irregular working conditions.

Any attempt to violate, adjust or repair this equipment by individuals not authorized by Instramed Indústria Médico Hospitalar Ltda. will automatically invalidate this warranty. This also applies in case of alterations made to this contract, the fiscal receipt, or to the equipment's serial number.

Instramed Indústria Médico Hospitalar Ltda. is not responsible for the improper use of this equipment, by people who are not familiar with its function or the techniques recommend for its proper use.

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ATTENTION: Instramed assumes no responsibility for any damage caused to individuals or property brought by failure to use this product in accordance with the information, recommendations and warnings presented in the user manual, alterations made in the device, attempts of repair not provided by authorized technical assistance centers, operation by unqualified personnel, use of defective device or use of accessories and parts not supplied by the manufacturer.

For information about warranty or technical assistance, please contact Instramed's technical support.

Copyright © 2016 Instramed. The Apolus, Instramed and its respective logos are trademarks of Instramed Indústria Médico Hospitalar Ltda. The internal software of this product is Instramed's intellectual property, being protected under international copyright laws. It is provided exclusively to be used with this present device, identified by the serial number, and may not be, in whole or in part, evaluated, recompiled or altered in any way.

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