

INSTRAMED

CPR Maestro

Efficiency in saving lives.


The accessory for CPR feedback
that ensures greater safety and effectiveness in emergency care.



 **CPR**
Maestro

CPR feedback device

The CPR Maestro is an accessory that assists and guides the application of CPR maneuvers, making the process much more efficient and reliable, significantly increasing the chances of reversing cardiac arrest. Its use enables a more effective and consistent performance of CPR, resulting not only in higher patient survival rates but also reducing the occurrence of severe sequelae due to lack of oxygenation.



How it works

The RCP Maestro is positioned on the patient's chest and connected to the defibrillator or cardioverter, providing real-time feedback on the quality of CPR. Through voice prompts and visual indicators, it helps the rescuer maintain the ideal frequency and depth of chest compressions. Its dedicated display shows this information directly on the device, and when used with cardioverters, the data is also displayed on the equipment's screen, ensuring efficient support during care.



The CPR Maestro is compatible with:

i.on

i.on
pro

CARDIOMAX

DUALMAX

Also compatible with Eazy line.

Instramed Indústria Médico-Hospitalar Ltda.

RS Unit
Avenida Severo Dullius, 1395 - Sala 403
ZIP 90200-310 - Porto Alegre - RS - Brazil
Tel.: +55 51 2313 3430

SC Unit
Rua Albatroz, 237 - Pedra Branca
ZIP 88137-290 - Palhoça - SC - Brazil
Tel.: +55 48 4042 4427

Folder RCP Maestro R3.1 2025 - EN

WWW.INSTRAMED.COM

INSTRAMED